

Culture conditions:

Overview

Primary human GBM cells are cultured in serum-free medium supplemented with EGF and FGF and can be grown as either floating neurospheres¹ or as an adherent monolayer as glioma neural stem (GNS)^{2,3} cells on a basement membrane of laminin or matrigel.

Culture medium

The cells have been tested using two serum free culture conditions with equivalent growth patterns.

1. **RHB-A**® cat#Y40001 <http://www.clontech.com>
2. **StemPro**® **NSC SFM** cat #A1050901 www.thermofisher.com (use as per manufacturer's instructions).

Growth factor supplementation

Epidermal growth factor (EGF) - final concentration (20ng/ml)

Fibroblast growth factor (FGF) - final concentration (10ng/ml)

Basement membrane attachment

For adherent growth purified laminin cat#L2020 <http://www.sigmaaldrich.com> or Corning® Matrigel® cat#354234 <http://www.corning.com> (Corning) are diluted 1:100 in cold phosphate buffered saline (PBS) and then filter sterilised. Tissue culture flasks are pre-coated and left to stand in the incubator for a minimum of three hours. Aspirate PBS and seed cells at the desired density.

Cell dissociation

ACCUTASE® cat# A6964-100ML <http://www.sigmaaldrich.com>

Each bottle contains 100 ml of 1x ACCUTASE enzyme in Dulbecco's PBS (0.2 g/L KCL, 0.2 g/L KH₂PO₄, 8 g/L NaCl, and 1.15 g/L Na₂HPO₄) containing 0.5 mM EDTA•4Na and 3 mg/L Phenol Red.

Use as per manufacturer's instructions.

Cell differentiation

To differentiate cells, remove EGF and FGF from the culture media and replace with 2% foetal bovine serum (FBS). Cells will differentiate within a few hours and adhere to non-matrigel/laminin coated tissue culture flasks and stain positive for neuronal and glial differentiation markers.

References

1. Reynolds BA, Weiss S. Generation of neurons and astrocytes from isolated cells of the adult mammalian central nervous system. *Science* 1992;255:1707-10.

2. Pollard SM, Yoshikawa K, Clarke ID, et al. Glioma stem cell lines expanded in adherent culture have tumor-specific phenotypes and are suitable for chemical and genetic screens. *Cell Stem Cell* 2009;4:568-80.
3. Day BW, Stringer BW, Wilson J, et al. Glioma surgical aspirate: a viable source of tumor tissue for experimental research. *Cancers* 2013;5:357-371.